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plugs are intended for gravity systems. Combi plug is not intended for administration of medicines.

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Warning :

		<p>Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.</p> <p>Period of use : The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the combi plug according to national guidelines.</p> <p>Disposal: Disposal according to local guidelines and/or clinical protocols.</p> <p>Note for users: Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.</p> <p>Storage: Products must be stored in the original manufacturer's packaging. Storage temperature must be +5°C to +30°C and relative humidity must be max. 75%.</p> <p>Important Notice: Date of last revision of the instructions for use: 9/5/2024</p>
LV	Latvia	<p>REF:</p> <p>CPHF-WYJ32R CPHF-WYJ32B CPHF-WYJ32W</p> <p>Svarīgs paziņojums: Lūdzu, izlasiet informāciju uz primārā iepakojuma. Nelietot, ja iepakojums ir bojāts. Sterilizēts ar EO.</p>

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CPHF-WYJ32R
CPHF-WYJ32B
CPHF-WYJ32W

Important notice :

Please read the information on the primary packaging. Do not use if package is damaged. Sterilized by EO.

Designated purpose:

Disposable medical devices designed for safety closure Luer closure. The combination of MLL and FLL connection connectors allows for safe and versatile use when closing syringes, intravenous catheters, three-way stopcocks, ramps and other infusion systems. Combi plugs are intended for gravity systems. Combi plug is not intended for administration of medicines.

Indications:

The combined MLL and FLL ensure safe closure of syringes, intravenous catheters, three-way stopcocks, ramps and other infusion systems. Combi plugs are intended for gravity systems.

Instructions for use:

1. while observing hygienic conditions, remove the stopper from the package.
2. Make the connection by twisting and making a secure connection.
3. Disconnect the plug by turning it counter-clockwise.

Patient population:

The product can be used by all patients. It has no gender or age restrictions.

Designated user:

The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).

Contraindications :

Medicines / solutions for which incompatibility with any of the materials used has been proven.

Residual risks / Side effects :

If the combi plug is used in accordance with the instructions given in the Intended Purpose section, no adverse reactions are known.

Safety measures :

Visual inspection of the Kombi stopper before use to ensure that the packaging is not damaged. Do not use if protective caps are missing or loose. Follow general safety precautions. Check the tightness of all connections at all times.

Warning :

Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.

Period of use :

The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the combi plug according to national guidelines.

Disposal:

Disposal according to local guidelines and/or clinical protocols.

Note for users:

Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.

Storage: Products must be stored in the original manufacturer's packaging. Storage temperature must be +5°C to +30°C and relative humidity must be max. 75%.

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