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2. Close the regulating button on the tube.
3. Remove the puncture tip cover and the LL connector cover.
4. Insert the puncture tip into the cap of the bottle or bag. Hang the bottle or bag on the rack.
5. Hand pump the drip chamber and fill it to 1/3, open the regulating clamp, open the vent valve and remove the air from the tubing.
6. Close the regulator.
7. Connect the set to the cannulas in the patient's body, open the pressure regulator and start the transfusion.

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		<p>Svarbus pranešimas: Nepradėkite perpylimo, jei pažeistas rinkinys arba kraujo maišelis. Visada laikykitės gamintojo nurodymų.</p> <p>Paskutinio naudojimo instrukcijos peržiūros data: 2024-09-05</p>
MT	Malta	<p>REF :</p> <p>CPHF40020 – Transfusion set Unohem DEHP FREE CPHF40023 – Transfusion set DuoHem DEHP FREE</p> <p>Important notice : Please read the information on the primary packaging. Do not use if package is damaged. Sterilized by EO.</p> <p>Designated purpose: Transfusion set for bottles and bags with blood. The single-use transfusion set is intended for administering blood to the circulatory system using intravenous cannulas.</p> <p>Indications: The single-use transfusion set is intended for administering blood to the circulatory system using intravenous cannulas.</p> <p>Instructions for use:            1. while observing hygienic conditions, remove the kit from the package.            2. Close the regulating button on the tube.            3. Remove the puncture tip cover and the LL connector cover.            4. Insert the puncture tip into the cap of the bottle or bag. Hang the bottle or bag on the rack.            5. Hand pump the drip chamber and fill it to 1/3, open the regulating clamp, open the vent valve and remove the air from the tubing.            6. Close the regulator.            7. Connect the set to the cannulas in the patient's body, open the pressure regulator and start the transfusion.</p> <p>Patient population:</p>

The product can be used by all patients. It has no gender or age restrictions.

**Designated user:**

The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).

**Contraindications :**

It must not be used other than as stated in the intended purpose. Do not use in patients with a known allergic reaction to any of the product's materials.

**Residual risks / Side effects :**

If the set is used in accordance with the instructions given in the Instructions for Use, no adverse reactions are known.

**Safety measures :**

Visual inspection of the infusion set before use for damage. – Do not use if protective caps are missing or loose.

**Warning :**

Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.

**Period of use :**

The duration of use depends on the treatment of the patient according to the summary of product information of the drug. Replace the application set according to national guidelines.

**Disposal:**

Disposal according to local guidelines and/or clinical protocols.

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