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**Important notice :**

Please read the information on the primary packaging. Do not use if package is damaged. Sterilized by EO.

**Designated purpose:**

The transfer needle is a double-sided transfer mandrel for mixing sterile solutions.

**Indications:**

dilution of lyophilisates, mixing of solutions and medicines

**Instructions for use:**

1. while observing hygienic conditions, remove the spike from the package.
2. Remove the protective covers from the mandrel.
3. Hold the lancing mandrel by the ergonomic wings, insert into the stopper of one bottle, and insert the other end of the mandrel into the stopper of the other bottle. Start mixing the solutions.

**Patient population:**

The product can be used in all patients who are prescribed infusion therapy. It has no gender or age restrictions.

**Designated user:**

The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).

**Contraindications :**

Do not use in patients with a known allergic reaction to any of the product's materials.

**Residual risks / Side effects :**

No side effects are known. The following complications may occur during the fluid transfer, dissolution, or dilution procedure

Microbiological or foreign particle contamination to infection due to insufficient aseptic handling.

Safety measures :

Do not sterilize. Always refer to the Summary of Product Characteristics for the mixing procedure.

Warning :

Do not reuse. Reusing a single-use device poses a potential risk to the patient or user. It can lead to contamination and/or deterioration of the product's functionality. Contamination or deterioration of the device's functionality can cause injury, illness, or death of the patient. Each transfer tip is intended for one use with one container and is not intended for reuse. Do not use the product if the packaging or the product is damaged or contaminated. The use of a damaged or contaminated product can cause injury, illness, or death of the patient. The duration of use is determined by the medicinal product used and its Summary of Product Characteristics and national guidelines or hospital protocols.

Period of use :

The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the transfer tip according to national guidelines.

Disposal:

Disposal according to local guidelines and/or clinical protocols.

Note for users:

Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.

Storage: Products must be stored in the original manufacturer's packaging. Storage temperature must be +5°C to +30°C and relative humidity must be max. 75%.

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3. Hold the lancing mandrel by the ergonomic wings, insert into the stopper of one bottle, and insert the other end of the mandrel into the stopper of the other bottle. Start mixing the solutions.

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**Residual risks / Side effects :**

No side effects are known. The following complications may occur during the fluid transfer, dissolution, or dilution procedure  
Microbiological or foreign particle contamination to infection due to insufficient aseptic handling.

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