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Instructions for use:

1. while observing hygienic conditions, remove the tube from the package.
2. Remove the protective covers from the connection connectors.
3. Connect the LLF connector of the tubing to the syringe and the LLM connector of the tubing to the catheter.

Patient population:

The product can be used in all patients who are prescribed infusion therapy. It has no gender or age restrictions.

Designated user:

The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).

Contraindications :

Medicines / solutions for which incompatibility with any of the materials used has been proven. The Summary of Product Characteristics of the respective drugs/solutions should be taken into account.

Residual risks / Side effects :

General risks of infusion therapy including air embolism, local or systemic infection, contamination by foreign particles. Another hazard to patients, visitors or staff may be the leakage of cytotoxic or biological fluids. Adverse reactions are not known if the tube is used in accordance with the instructions given in the Precautions section.

Safety measures :

Visual inspection of the perfusion tubing before use for damage. – Do not use if protective caps are missing or loose. Follow general safety precautions. Follow directions for use for parallel infusions. Please also take into account the information in the summary of product characteristics of the respective drugs/solutions regarding possible incompatibilities of these drugs/solutions. Check the tightness of all connections throughout the infusion.

Warning :

		<p>Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.</p> <p>Period of use : The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the perfusion tubing according to national guidelines.</p> <p>Disposal: Disposal according to local guidelines and/or clinical protocols.</p> <p>Note for users: Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.</p> <p>Storage: Products must be stored in the original manufacturer's packaging. Storage temperature must be +5°C to +30°C and relative humidity must be max. 75%.</p> <p>Important Notice: Date of last revision of the instructions for use: 9/5/2024</p>
LV	Latvia	<p>REF:</p> <p>CPHF50150 / CPHF52025 / CPHF52250 CPHF50200 / CPHF52100 / CPHF52300 CPHF50300 / CPHF52150 CPHF52010 / CPHF52200</p> <p>Svarīgs paziņojums:</p>

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Warning :

		<p>Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.</p> <p>Period of use : The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the perfusion tubing according to national guidelines.</p> <p>Disposal: Disposal according to local guidelines and/or clinical protocols.</p> <p>Note for users: Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.</p> <p>Storage: Products must be stored in the original manufacturer's packaging. Storage temperature must be +5°C to +30°C and relative humidity must be max. 75%.</p> <p>Important Notice: Date of last revision of the instructions for use: 9/5/2024</p>
PL	Poland	<p>ODN.:</p> <p>CPHF50150 / CPHF52025 / CPHF52250 CPHF50200 / CPHF52100 / CPHF52300 CPHF50300 / CPHF52150 CPHF52010 / CPHF52200</p> <p>Ważna uwaga:</p>

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