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Important notice :

Please read the information on the primary packaging. Do not use if package is damaged. Sterilized by EO.

Designated purpose:

Single-use sterile infusion sets for infusion/drug administration - marked with the symbol "G" (gravity, gravity) on the primary package. Single-use sterile infusion systems for the administration of infusion/drug by fall administration and compatible pumps (2 bar pressure resistant) - marked with the symbol "P" (pressure) on the primary packaging.

Indications:

Infusion sets are used for infusion therapy, e.g.: maintenance of fluid volume, fluid exchange, infusion of total parenteral nutrition, infusion of drugs (e.g. antibiotics) using bottles and bags, and administration of chemotherapy drugs.

Instructions for use:

1. while observing hygienic conditions, remove the set from the packaging.
2. Close the regulator on the tubing.
3. Remove the puncture tip cover and the LL connector cover.
4. Insert the puncture tip into the cap of the bottle or bag. Hang the bottle or bag on the rack.
5. Pump the drip chamber by hand and fill it to 1/3, open the regulator, open the vent valve and remove the air from the tubing.
6. Close the regulator.
7. Connect the set to the intravenous cannula in the patient's body, open the control clamp and start the infusion.

Patient population:

The product can be used in all patients who are prescribed infusion therapy. It has no gender or age restrictions.

Designated user:

The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).

	<p>Contraindications: Medicines / solutions for which incompatibility with any of the materials used has been proven. Blood products or blood components must not be administered. The Summary of Product Characteristics of the respective drugs/solutions should be taken into account.</p> <p>Residual risks / Side effects : General risks of infusion therapy including air embolism, local or systemic infection, contamination by foreign particles. Another hazard to patients, visitors or staff may be the leakage of cytotoxic or biological fluids.</p> <p>Safety measures : Visual inspection of the infusion set before use for damage. – Do not use if protective caps are missing or loose. When using a positive pressure infusion, first remove the air from the container. Follow general safety precautions. Follow directions for use for parallel infusions. Please also take into account the information in the summary of product characteristics of the respective drugs/solutions regarding possible incompatibilities of these drugs/solutions. Check the tightness of all connections throughout the infusion. – Always follow the instructions for use of the manufacturer of infusion pumps (overpressure infusion sets). The instructions for use of the infusion pump manufacturer must be followed.</p> <p>Warning : Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.</p> <p>Period of use : The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the application set according to national guidelines. The duration of use of intravenous administration sets that can be used for pressure infusion is limited to 24 h when used with an infusion pump.</p> <p>Disposal: Disposal according to local guidelines and/or clinical protocols. For disposal, use the integrated tip protection on the control terminal.</p>
--	---

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	<p>Single-use sterile infusion sets for infusion/drug administration - marked with the symbol "G" (gravity, gravity) on the primary package. Single-use sterile infusion systems for the administration of infusion/drug by fall administration and compatible pumps (2 bar pressure resistant) - marked with the symbol "P" (pressure) on the primary packaging.</p> <p>Indications: Infusion sets are used for infusion therapy, e.g.: maintenance of fluid volume, fluid exchange, infusion of total parenteral nutrition, infusion of drugs (e.g. antibiotics) using bottles and bags, and administration of chemotherapy drugs.</p> <p>Instructions for use:</p> <ol style="list-style-type: none"> 1. while observing hygienic conditions, remove the set from the packaging. 2. Close the regulator on the tubing. 3. Remove the puncture tip cover and the LL connector cover. 4. Insert the puncture tip into the cap of the bottle or bag. Hang the bottle or bag on the rack. 5. Pump the drip chamber by hand and fill it to 1/3, open the regulator, open the vent valve and remove the air from the tubing. 6. Close the regulator. 7. Connect the set to the intravenous cannula in the patient's body, open the control clamp and start the infusion. <p>Patient population: The product can be used in all patients who are prescribed infusion therapy. It has no gender or age restrictions.</p> <p>Designated user: The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).</p> <p>Contraindications: Medicines / solutions for which incompatibility with any of the materials used has been proven. Blood products or blood components must not be administered. The Summary of Product Characteristics of the respective drugs/solutions should be taken into account.</p> <p>Residual risks / Side effects :</p>
--	---

General risks of infusion therapy including air embolism, local or systemic infection, contamination by foreign particles. Another hazard to patients, visitors or staff may be the leakage of cytotoxic or biological fluids.

Safety measures :

Visual inspection of the infusion set before use for damage. – Do not use if protective caps are missing or loose. When using a positive pressure infusion, first remove the air from the container. Follow general safety precautions. Follow directions for use for parallel infusions. Please also take into account the information in the summary of product characteristics of the respective drugs/solutions regarding possible incompatibilities of these drugs/solutions. Check the tightness of all connections throughout the infusion. – Always follow the instructions for use of the manufacturer of infusion pumps (overpressure infusion sets). The instructions for use of the infusion pump manufacturer must be followed.

Warning :

Reuse of single-use devices poses a potential risk to the patient or user. It can lead to contamination or deterioration of the device's functionality, which can cause injury, illness, or death of the patient. Administration of light-sensitive drugs requires the use of UV-protected intravenous administration sets.

Period of use :

The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the application set according to national guidelines. The duration of use of intravenous administration sets that can be used for pressure infusion is limited to 24 h when used with an infusion pump.

Disposal:

Disposal according to local guidelines and/or clinical protocols. For disposal, use the integrated tip protection on the control terminal.

Note for users:

Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.

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