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Designated purpose:

Coro-Spike is a vented collection cannula (aspiration spike) for collecting fluids from glass/plastic vials.

Indications:

Collection and injection from and into multi-dose vials, and subsequent application to the infusion line

Preparation of cytostatic injections

Operating instructions

1) While maintaining hygienic conditions, remove the protective cover from the medicine container and disinfect the cap.

2) Remove the protective cover from the tip

3) Insert the tip.

4) Open the cover.

5) Connect the syringe and aspirate/draw the desired volume.

6) Important: The solution from the syringe into the container must be injected upside down. Risk of excessive pressure and impaired venting.

Patient population:

The product can be used in all patients who are prescribed infusion therapy. It has no gender or age restrictions.

Designated user:

The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).

Contraindications :

Collection of cellular material containing blood samples. Dissolving solid medicinal products with another solvent. Special transmission means are available for this indication

Residual risks / Side effects :

No side effects are known. The following complications may occur during the fluid transfer, dissolution, or dilution procedure

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	<p>Designated purpose: Coro-Spike is a vented collection cannula (aspiration spike) for collecting fluids from glass/plastic vials.</p> <p>Indications: Collection and injection from and into multi-dose vials, and subsequent application to the infusion line Preparation of cytostatic injections</p> <p>Operating instructions 1) While maintaining hygienic conditions, remove the protective cover from the medicine container and disinfect the cap. 2) Remove the protective cover from the tip 3) Insert the tip. 4) Open the cover. 5) Connect the syringe and aspirate/draw the desired volume. 6) Important: The solution from the syringe into the container must be injected upside down. Risk of excessive pressure and impaired venting.</p> <p>Patient population: The product can be used in all patients who are prescribed infusion therapy. It has no gender or age restrictions.</p> <p>Designated user: The device is intended for use by authorized persons, eg: health care professionals, licensed and/or board-certified physicians, medical professionals, and/or specially trained patients or caregivers (according to local regulations).</p> <p>Contraindications : Collection of cellular material containing blood samples. Dissolving solid medicinal products with another solvent. Special transmission means are available for this indication</p> <p>Residual risks / Side effects : No side effects are known. The following complications may occur during the fluid transfer, dissolution, or dilution procedure</p>
--	---

Microbiological or foreign particle contamination to infection due to insufficient aseptic handling.

Safety measures :

Do not sterilize. Always look up the mixing procedure in the Product Summary.

Warning :

Do not reuse. Reusing a single-use device poses a potential risk to the patient or user. It can lead to contamination and/or deterioration of the product's functionality. Contamination or deterioration of the device's functionality can cause injury, illness, or death of the patient. Each dispensing tip is intended for one use with one container and is not intended for reuse. Do not use the product if the packaging or the product is damaged or contaminated. The use of a damaged or contaminated product can cause injury, illness, or death of the patient. The duration of use is determined by the medicinal product used and its Summary of Product Characteristics and national guidelines or hospital protocols.

Period of use :

The duration of use depends on the treatment of the patient according to the Summary of Product Characteristics of the drug/solution. Replace the application set according to national guidelines.

Disposal:

Disposal according to local guidelines and/or clinical protocols.

Note for users:

Any serious incident that occurs in connection with this device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.

Storage: Products must be stored in the original manufacturer's packaging. Storage temperature must be +5°C to +30°C and relative humidity must be max. 75%.

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